



BiOfactor MTA

Root Canal Filling and Repair Material




Easy application

No special equipment required for mixing or placement.

Fast setting

Finer powder for faster hydration.






What is Biofactor MTA ?


BiofactorMTA(mineral trioxide aggregate) is a tricalcium silicate-based white coloured bioactive repair cement that can be used universally endodontic and pediatric indications.

- BiofactorMTA kits contain gel liquid that improves the handling and placement.
- The gel liquid enables you to use the material with different consistencies (thick, putty and flow) for a variety of procedures.
- Finer powder for faster hydration
- Better mixing, handling and placement(based on procedure flow or puty consistency)
- High Radiopacity
- No special equipment required for mixing or placement
- Biocompatible
- Strong Sealing Properties
- No Bismuth,won't discolor teeth
- Thin consistency can be precisely delivered into the canal with a tips for apexification and apical plug.




Flow Consistency

All of the powder is mixed with 3 drops of liquid. The syringe is filled with a spatula and applied.



Putty Consistency

All of the powder is mixed with 2 drops of liquid. The mixed cement applied with the a spatula.



Thick Consistency

Take a scoope of powder from brown glass vial with using big direction of spoon, mixed with 4 drops of liquid. The syringe is filled with a spatula and applied.



RESOPTION
REZOPSIYON

APEXIFICATION
APEKSİFİKASYON

APICAL PLUG
APIKAL TIKAÇ



PULP CAPPING
PULPA KUAFAJI

PULP PERFORATION
PULPA PERFORASYONU

PULPOTOMY
PULPATOMİ

ROOT END FILLING
KÖK BİTİM DOLGUSU

Indications

Flow Consistency

- Resorption
- Apexification
- Apical plug


Indications

Putty and Flow Consistency

- Pulp capping
- Pulp chamber perforation
- Pulpotomy
- Root end filling.



1.5-year follow-up after repair of internal resorption and apical lesion with BioFactor MTA.




1.5-year follow-up after repair of a traumatized open apex lesion with BioFactor MTA.



6 months follow-up after repair of traumatized permanent teeth with BioFactor MTA.


Shear Bond Strength Of Two Calcium Silicate-Based Cements To Compomer

Burak BULDUR, Fatih ÖZNURHAN, Mevlut KAYABAŞI, Feride ŞAHİN
Cumhuriyet Dental Journal: 2018; 21(1)



The Push-Out Bond Strength of BIOfactor Mineral Trioxide Aggregate, A Novel Root Repair Material

Makbule Bilge AKBULUT, Durmuş Alperen BOZKURT, Arslan TERLEMEZ, Melek AKMAN
The Korean Academy of Conservative Dentistry: 2019



No statistically significant differences in the push-out bond strength were found between the three groups of materials ($p > 0.017$).

